




PULLUP®

The Versatile Adjustable Loop
Suspensory Fixation Device

Visit us at
Stand 22



Including:



Pull, flip, adjust

Solution for cortical fixation of knee ligament reconstructions:

Reproducible


Controlled locking at the end of the blind tunnel - bone preserving

Simple


No preoperative planning of tunnel length - intraoperative adjustment

Reliable


Stable locking system with renowned materials



A system that can be adapted to your surgical approach



For cortical tunnels with \varnothing 4.5 mm



For cortical tunnels with \varnothing from 5 to 10 mm

Speak to your local Sales Specialist for further information or contact us using the details below:

T: 01443 719 555

E: info@hospitalinnovations.co.uk

www.hospitalinnovations.com

Hospital Innovations Limited
Concept House
Talbot Green Business Park
Pontyclun
CF72 9FG

HOSPITAL
INNOVATIONS