

# ReadiGRAFT® Cancellous Chips



**ReadiGRAFT® Cancellous Chips** are a natural osteoconductive scaffold. Chips are available in 4 grind sizes and a variety of volumes to provide maximum surgical flexibility.

<p>Clinical Applications</p>	<p>Any surgical application that requires bone void filler</p>
<p>Features And Benefits</p>	<p><b>Natural Osteoconductive Scaffold:</b> The natural interconnected porous architecture provides an optimal scaffold for cellular attachment, proliferation, and vascular in-growth.</p> <p><b>Hydrophilic:</b> Absorbs and retains bioactive fluids like blood, Platelet Rich Plasma (PRP), and Bone Marrow Aspirate (BMA).</p> <p><b>Sterile:</b> The implant is sterilised using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level of 10<sup>-6</sup>, without compromising the implant's inherent osteoconductivity or osteoinductive potential.</p>
<p>Alternative Products</p>	<p>Cancellous Cubes</p> <p>Cortical Cancellous Chips</p> <p>Femoral Head (Grinder)</p>

# ReadiGRAFT®

## Cancellous Chips

### CANCELLOUS CHIPS

Grind Size	Volume	Preservon Order No.	Freeze Dried Order No.	Frozen Order No.
0.1-2 mm	5 cc	PCAN5 12		
0.1-2 mm	15 cc	PCAN15 12		
0.1-2 mm	30 cc	PCAN30 12		
1-4 mm	5 cc	PCAN5 14		
1-4 mm	15 cc	PCAN15 14	CAN15 14BP	FCAN15 14
1-4 mm	30 cc	PCAN30 14	CAN30 14BP	FCAN30 14
1-4 mm	60 cc	PCAN60 14	CAN60 14BP	
1-4 mm	90 cc	PCAN90 14	CAN90 14BP	
1-8 mm	5 cc	PCAN5	CAN5	
1-8 mm	15 cc	PCAN15	CAN15	FCAN15
1-8 mm	20 cc	PCAN1/4	CAN1/4	
1-8 mm	30 cc	PCAN30	CAN30	FCAN30
1-8 mm	40 cc	PCAN1/2	CAN1/2	FCAN40
1-8 mm	60 cc	PCAN60	CAN60	
1-8 mm	80 cc			FCAN80
1-8 mm	90 cc	PCAN90	CAN90	FCAN9
4-10 mm	15 cc	PCAN15 410	CAN15 410BP	FCAN15 410
4-10 mm	30 cc	PCAN30 410	CAN30 410BP	FCAN30 410
4-10 mm	60 cc	PCAN60 410	CAN60 410BP	FCAN60 410
4-10 mm	90 cc	PCAN90 410	CAN90 410BP	FCAN90 410

Speak to your local Sales Specialist for further information or contact us using the details below:

T: 01443 719 555  
E: info@hospitalinnovations.co.uk

[www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited  
Concept House  
Talbot Green Business Park  
Pontyclun  
CF72 9FG

HOSPITAL  
INNOVATIONS