



# PULLUP®

The **Versatile** Adjustable Loop  
Suspensory Fixation Device




**Inside-Out**  
PULLUP® & PULLUP® XL




**All Inside**  
PULLUP® on both

Including:



**GraftTech®**  
Preparation station



**Outside-In**  
PULLUP® XL on both

# Pull, flip, adjust

Solution for cortical fixation of knee ligament reconstructions:

## Reproducible


Controlled locking at the end of the blind tunnel - bone preserving

## Simple


No preoperative planning of tunnel length - intraoperative adjustment

## Reliable


Stable locking system with renowned materials



A system that can be adapted to your surgical approach



For cortical tunnels with  $\varnothing$  4.5 mm



For cortical tunnels with  $\varnothing$  from 5 to 10 mm

Speak to your local Business Development Manager for further information or contact us using the details below:

T: 01443 719 555

E: [customerservice@hiuk.co.uk](mailto:customerservice@hiuk.co.uk)

[www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited

Concept House

Talbot Green Business Park

Pontyclun

CF72 9FG

