

ReadiGRAFT®

Cancellous Chips



ReadiGRAFT® Cancellous Chips are a natural osteoconductive scaffold. Chips are available in 4 grind sizes and a variety of volumes to provide maximum surgical flexibility.

Clinical Applications	Any surgical application that requires bone void filler
Features And Benefits	<p>Natural Osteoconductive Scaffold: The natural interconnected porous architecture provides an optimal scaffold for cellular attachment, proliferation, and vascular in-growth.</p> <p>Hydrophilic: Absorbs and retains bioactive fluids like blood, Platelet Rich Plasma (PRP), and Bone Marrow Aspirate (BMA).</p> <p>Sterile: The implant is sterilised using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level of 10⁻⁶, without compromising the implant's inherent osteoconductivity or osteoinductive potential.</p>
Alternative Products	<p>Cancellous Cubes</p> <p>Cortical Cancellous Chips</p> <p>Femoral Head (Grinder)</p>

ReadiGRAFT®

Cancellous Chips

Freeze-dried Chips

Item No.	Item Description
CAN15 14BP	Freeze-dried Crushed Cancellous Chips: 1-4mm, 15cc
CAN30 14BP	Freeze-dried Crushed Cancellous Chips: 1-4mm, 30cc
CAN60 14BP	Freeze-dried Crushed Cancellous Chips: 1-4mm, 60cc
CAN5	Freeze-dried Crushed Cancellous Chips: 1-8mm, 5cc
CAN15	Freeze-dried Crushed Cancellous Chips: 1-8mm, 15cc
CAN30	Freeze-dried Crushed Cancellous Chips: 1-8mm, 30cc
CAN60	Freeze-dried Crushed Cancellous Chips: 1-8mm, 60cc
CAN15 410BP	Freeze-dried Crushed Cancellous Chips: 4-10mm, 15cc
CAN30 410BP	Freeze-dried Crushed Cancellous Chips: 4-10mm, 30cc
CAN60 410BP	Freeze-dried Crushed Cancellous Chips: 4-10mm, 60cc
CAN90 410BP	Freeze-dried Crushed Cancellous Chips: 4-10mm, 90cc

Frozen Chips

Item No.	Item Description
FCAN30 14	Frozen Crushed Cancellous Chips: 1-4mm, 30cc
FCAN60 14	Frozen Crushed Cancellous Chips: 1-4mm, 60cc
FCAN30	Frozen Crushed Cancellous Chips: 1-8mm, 30cc

Speak to your local Business Development Manager for further information or contact us using the details below:

T: 01443 719 555

E: customerservice@hiuk.co.uk

www.hospitalinnovations.com

Hospital Innovations Limited
Concept House
Talbot Green Business Park
Pontyclun
CF72 9FG

