



OraGraft®

Cancellous Mineralised Powder

Clinical Overview The OraGraft line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

- Applications**
- Periodontal Defects
 - Ridge Maintenance
 - Extraction Site Preservation
 - Sinus Grafting
 - Craniofacial Reconstruction

- Features & Benefits**
- **Increased Safety:** Sterile with an SAL 10^{-6} .
 - **Convenient Storage:** Ambient temperature with three-year shelf life.
 - **Graft Incorporation:** Osteoconductive potential with a combination of open trabecular architecture and dense particles allows for healing through “creeping substitution.”

Speak to your local Business Development Manager for further information or contact us using the details below:

T: 01443 719 555

E: customerservice@hiuk.co.uk

www.hospitalinnovations.com

Hospital Innovations Limited

Concept House

Talbot Green Business Park

Pontyclun

CF72 9FG





OraGraft Cancellous Mineralised Powder

Ambient Temperature*/3 year shelf life

250 - 1000 microns

Size	Order Code
0.50 cc	OCAN-0.5A
1.00 cc	OCAN-1.0A
2.00 cc	OCAN-2.0A

1000 - 2000 microns

Size	Order Code
0.50 cc	OCAN-0.5B
1.00 cc	OCAN-1.0B
2.00 cc	OCAN-2.0B

*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)

Speak to your local Business Development Manager for further information or contact us using the details below:

T: 01443 719 555

E: customerservice@hiuk.co.uk

www.hospitalinnovations.com

Hospital Innovations Limited

Concept House

Talbot Green Business Park

Pontyclun

CF72 9FG

