


# NanoFx®

## A-Curve

### Smaller. Deeper. Better.

"The overriding finding of the present study is that small subchondral drill holes reflecting the physiological subchondral trabecular distance significantly improve osteochondral repair."\*

*Eldracher M, et. al. Small subchondral drill holes improve marrow stimulation of articular cartilage defects. AJSM 2014 Nov;42(11):2741-50.*



- Specially designed curve for hard to access ankle lesions
- Reduces surface skyving for off-axis channel preparation
- Improves bone marrow access for cartilage repair procedures



Thumble



Handle Instrument




PleuriStik™ Guide Wire

Item No.	Item Description
FURS-010R	NanoFx® Thumb Tab Accessory
5500-4010	NanoFx® A Curve Handle
5500-1020	NanoFx® Handle 15
FURS-2101	NanoFx® PleuriStik 15 - Box of 5
FURS-4101	NanoFx® PleuriStik A - Box of 5
FURS-2101/1	NanoFx® PleuriStik 15 - Single
FURS-4101/1	NanoFx® PleuriStik - Single


## Smaller

"Significant enhancements were observed at the level of individual parameters and of overall histological articular cartilage repair, together with improved immunoreactivity to type II and type I collagen of the cartilaginous repair tissue. Second, the microarchitecture of both the subchondral bone plate and the subarticular spongiosa was better reconstituted."

*Eldracher M, Orth P, Cucchiari M, Pape D, Madry H. Small subchondral drill holes improve marrow stimulation of articular cartilage defects. Am J Sports Med. 2014 Nov;42(11):2741-50.*




## DEEPER



STEP 1

Assemble & orient on lesion



STEP 2


Gently tap to depth stop



STEP 3

Push thumb to extract pleuristik.  
Repeat.

## BETTER



Figures: open trabecular channels;  
 closed trabecular channels, microCT  
comparison: Axial (top), Sagittal (bottom).

Speak to your local Business Development Manager, for further information using the details below:

**T:** 01443 719 555

**E:** [customerservice@hiuk.co.uk](mailto:customerservice@hiuk.co.uk)

**W:** [www.hospitalinnovations.com](http://www.hospitalinnovations.com)

Hospital Innovations Limited  
Concept House  
Talbot Green Business Park  
Pontyclun  
CF72 9FG